

Connecting a hearing aid to a computer/laptop


It is possible to connect hearing aids directly to the audio out/headphone socket on a computer/laptop/tablet

You will need 'direct input leads'

These can be purchased here

<https://www.connevans.co.uk/product/1852/DPV600/600mm-blue-plugged-personal-stereo-V-lead-for-2-hearing-aids> for Phonak hearing aids

<https://www.connevans.co.uk/product/410/FMGV600/600mm-fmGenie---Oticon-binaural-V-direct-input-lead> for Oticon hearing aids



You will also need *audioshoes*, which are specific to your model of hearing aid(s). Use this link below to help you identify which audioshoe(s) you will need

<https://www.connevans.co.uk/catalogue/10888618/My-Hearing-Aid>

follow the steps until you get to *Direct Audioshoe Type*- it will be AS followed by a number eg AS18

This video shows you how to put the audioshoe onto the hearing aid

<https://www.youtube.com/watch?v=Yo4NmxGudL4>

(instead of adding the small receiver at the end, plug in the direct input leads)

Then plug into the headphone jack on your device



If you have a newer model of iPhone/iPad, you will need a lightning to 3.5mm headphone jack adapter



<https://www.apple.com/uk/shop/product/MMX62ZM/A/lightning-to-35mm-headphone-jack-adapter>